



**CONTIPHARMA**

*a member of Bedelco group*

**PROTECT**

## Masks FFP2 / KN95



OUTSIDE



INSIDE

**BRAND**

PURISM

**BACTERIAL FILTRATION EFFICIENCY**

≥ 95%



## SPECIFICATIONS

5-layer mask :

- First layer: 60g/m<sup>2</sup> of unwoven fabric
- Second and third layer: 30g/m<sup>2</sup> e high-efficiency unwoven fabric, blown by fusion
- Fourth layer: 50g/m<sup>2</sup> hot air cotton
- Fifth layer: 25g/m<sup>2</sup> of unwoven fabric

Composition :

Non-woven wire-bound, non-woven smelted, hot air cotton, nose clip, earrings

Product performance :

perfect fit for each face and optimal protection. Suitable for the protection of certain non-oil particles, such as flu virus, dust, pollen, etc.

Recommended port time :

4-8 hours, to be changed between medical procedures. After contact with the used mask, be sure to disinfect your hands.

---

## APPLICATION

Schools, high-speed trains, airports, public transport and other relatively densely populated places and family environment.

---

## COMPLIANCE

EC marking according to European Directive 89/686/EEC

Standard 149:2001 - A1:2009 / GB2626-2006 KN95 in connection with REGULATION CE R2016/425 (Personal Protective Equipment)



## CONDITIONING

5 pcs/bag, 4 bags/box, 48 boxes/cardboard, 960 pieces/cardboard,  
 Cardboard size: 57.5 × 33.5 × 38 cm - 8.35 kg

## STORAGE

Store in a cool, dry and clean place, away from the fire, between -20 and 38 degrees C with a relative humidity of less than 30%.

## STORAGE TIME

3 years