

Nitrile gloves « Purism » G001

Ref 50001



SPECIFICATIONS

G001 «Purism» examination gloves are suitable for work environments where there is possible contact with bodily fluids, microorganisms and chemicals. They do not contain natural rubber latex and are an excellent alternative for people with Type I allergies.

These are non-powder gloves that offer great flexibility, high sensitivity and resistance to chemicals. They are tear-resistant, ambidextrous and non-sterile.

SIZE				
Size	Minimal	Palm	Weight	Minimum
	length (mm)	width (mm)	(g)	thickness (mm)
XS	240	≤80	3.8 ± 0.2	0,05
S	240	80 ± 10	4.3 ± 0.2	0,05
Μ	240	95 ± 10	4.8 ± 0.2	0,05
L	240	110 ± 10	5.3 ± 0.2	0,05
XL	240	≥110	5.8 ± 0.2	0,05

STANDARD

EN455





MAIN COMMODITIES

- · Butadiene nitrile rubber ≥90%
- · Zinc oxide 4%
- · Titanium dioxide 2.6%

APPLICATION

Medical, food, chemical industries,...

REMARKS

- · Examination gloves are single-use products and should only be used for one procedure and/or one patient
- \cdot When wearing gloves, do not wear rings or other ornaments and cut nails gently
- · Wear them indifferently on the left or right hand, please choose gloves adapted to the specifications of your hand
- $\cdot\,$ Do not use if gloves are contaminated due to defective packaging
- · Do not directly irradiate strong light such as sunlight or ultraviolet rays.

CONDITIONING

100-piece cardboard.

STORAGE

Gloves should be stored in a cool, dry and dark environment without corrosive gases.

Belgium Web : www.contipharma-protect.com